

Cleaning and Maintenance

A Note about the Upholstery

All Knight chair upholstery is made from the highest quality materials available and is designed for maximum patient comfort, easy cleaning and durability. Knight chair upholstery will provide years of service when correctly maintained according to the Operator and Care Manual. Like all fine upholstery some minor wrinkles, waves or changes may occur. These changes are natural and in no way effect the performance or life of the product.

Cleaning Instructions - Upholstery

When cleaning upholstery, it is suggested that you test all recommended commercial cleaning agents in an inconspicuous spot to check for potential damage to the upholstery material before using on the entire area to be cleaned. For everyday cleaning we recommend a mild dishwashing liquid, warm water and a damp cloth. Spraying upholstery with a silicone-based coating will improve luster and soil resistance.

Cleaning Instructions - Painted Surfaces

All metal surfaces are powder-coated and resistant to scratching and scuffing. Clean with a clean cloth dampened with mild soapy water.

Recommended Infection Control Products

The disinfectant/cleaner products that we have tested and recommend for our upholstery and painted items are as follows:

1. CAVACIDE - distributed by E & D Dental Products, Inc., Westfield, NJ
2. PRECISE - Hospital Foam Cleaner Disinfectant - manufactured by Caltech Industries, Inc., Midland, MI

CAUTION: *Follow manufacturer's directions for concentration and application of disinfecting material. Avoid prolonged application of any disinfectant/cleaner products, because they may cause staining of upholstery material.*

Knight®
by Midmark

Midmark Corporation • P O Box 286 • 60 Vista Drive
Versailles • Ohio • 45380-0286 • www.midmark.com
Phone: 937-526-3662 Fax: 937-526-5542

TP1142 Rev. C